

UGANDA NATIONAL TUBERCUSLOSIS REFERENCE LABORATORY GENEXPERT MTB/RIF PT
SCHEME - DRIED TUBE SPECIMENS (DTS) Xpert® MTB/RIF ONLINE REPORTING
INSTRUCTIONS


PT MANAGEMENT SYSTEM GUIDE (Website: <https://srlugpt.com/>)

User NAME: «User» Password: «Password»

DO not share this username and password with other laboratories. Welcome to NTRL PT management system.


Access.

- i. The system can be accessed using a computer that is connected to internet.
- ii. Open browser (Chrome and firefox recommended) and enter the following address 'srlugpt.com' in your address bar and press Enter



Login.

From the home page click 'Login' from the top menu and the login page will load.
Enter your credentials (username and password as highlighted above)




Home page.

After login the system will take you to the home page with the profile of your facility



- 1 Account creation.
- 2 Scheme enrollment.




Account Profile

Result submission

To enter result, go to **Results** from top menu and choose 'Genexpert' this will load a panels' page

List of schemes



SCHEME - DRIED TUBE SPECIMENS (DTS) Xpert® MTB/RIF ONLINE REPORTING
INSTRUCTIONS

From the panels' page choose the **round** for which results you want to enter. This will load a list of panels **page**

List of Panels under Genexpert

GeneXPert PT ROUND TWO
2021-GX
2021 round 2

Choose Add PT results' from the sub menu (**Please first scan the result sheet and save it either as pdf or image**)

1 Results.

View entered resultsView result sheetAdd pt results

Add Genexpert results for GeneXPert%20PT%20ROUND%20TWO%202021

Assay model*

Enter your results and click on 'add results' button.

When you are done entering results for all the five panels, attach the results file and submit

1 Results.

View entered resultsView result sheetAdd pt results

GeneXpert results

10search

Country	Lab No.	Lab Name	PT-ID No.	Gx Serial No.	Module Model.	Module No.	Catr exp date.	Panel Received Date	Results Report Date	Date Tested	Panel Id/	MTB Detected	Rif Resistan
---------	---------	----------	-----------	---------------	---------------	------------	----------------	---------------------	---------------------	-------------	-----------	--------------	--------------

Attach Genexpert results for GeneXPert%20PT%20ROUND%20TWO%202021

Attached file

Browse


Add results


To view your entered results, click ‘**View entered Results**’


To view your performance report.

Choose ‘**Reports**’ from the top menu. **Note** the report will only be accessible after it’s been authorized by PT provider this may take some time from the time of submission.

Your Performance Report

 Feedback report pg 1 .

 Feedback report pg 2.

 Summary report.

Final Individual Participant Report

GeneXpert Proficiency Testing Scheme

GeneXPert%20PT%20ROUND%20TWO%202021

Date PT Shipped: **Sept. 6, 2021** Date PT received: **Sept. 7, 2021** Closing Date: **Oct. 18, 2021** Results return to PT provider: **Feb. 23, 2022** Results TAT 169 days, Report Issue date: Feb. 25, 2022

Country/Region

Site Name

Participant Code

The report can be printed and filed

PT 011 F5 Authorized by: PT manager